

EPI-SPORT® EPICONDYLITIS CLASP

Model No: 19-600

Sizes: XSmall - XLarge

Color: Black



INDICATIONS:

- Epicondylitis
- Tennis Elbow
- Muscle Strains
- Carpal Tunnel Syndrome
- Arthritis
- Tendinitis

BENEFITS:


This clasp helps reduce pain and discomfort on both sides of the forearm. The clasp applies dual compression on the long tendons of the forearm extensor and flexor muscles making them more parallel to allow healing to begin and reduce inflammation and pain.

PRODUCT FEATURES:

- Comfortable shock absorbing and moisture wicking liner
- Dual action compression on lateral and medial sides of elbow
- Elastic tension strap holds support in place

HOW TO APPLY:

Position on arm with strap hanging. Slip clasp on arm and position the large side on top with buckle facing inward. Position so large portion is directly over the extensor tendon. Tighten and secure strap until clasp is snug and comfortable. Readjust if necessary.



HOW TO APPLY:

With arm band loose and Velcro® strap closed, slip EpiSport® on arm and position the large side of the clasp on top with the buckle facing inward. Position EpiSport® where compression is desired and readjust strap for a snug fit. See enclosed instructions.




CARE INSTRUCTIONS:

Fasten closures. Hand wash in cold water with mild soap. Wrap in a towel to remove excess water (do not twist or wring). Air dry away from heat. Do not tumble dry. Do not bleach.

SIZING:

Measure forearm 1 1/2" to 2" below elbow.

Size	Fits
X-Small	8" - 9"
Small	9" - 10"
Medium	10" - 11"
Large	11" - 12"
X-Large	12" - 13 1/2"



800-327-4110